



EUROStar II Plus – a well-priced fluorescence microscope. Now available with LED and transmitted-light illuminator!



Technical Data

Dimensions

W x D x H approx. 200 x 350 x 450 mm
Weight approx. 7.2 kg

Power supply

Output voltage 12V (Battery use possible)
Supply voltage 100 to 240V
Power supply range autom. voltage conversion

LED-Light source EUROStar-Bluelight

LED type Azul EU1
Excitation light source wave length 460-490 nm
LED voltage 3.75-3.99V
Power 5W
Constant light flux at 460-490 nm 30 lm
Product life approx. 15,000 h
Laser classification 2M
Deterioration indicator acoustic alarm

LED transmitted-light

brightfield illuminator Carl Zeiss/EUROIMMUN

Filter set for FITC

Excitation filter/Emission filter 450-490 nm/515 nm
Beam splitter 510 nm

Opto-mechanical components

Objective change revolver manual, 4-fold
Objective 1 A-Plan 20x/0.45
Objective 2 A-Plan 40x/0.65
Eyepiece PL 10x/20 Br. and PL 10x/20 Br. foc.
Binocular phototube 30°/20
Sliding prism 100% vis / 100% doc
Max. number of viewing fields 20
Interpupillary distance 55 to 75 mm
Viewing angle/height 30°/430 mm
Visual port Tube factor 1x
Photo/Video port Tube factor 1x

Stand

Stage focusing coarse drive 4 mm/U
Stage focusing fine drive 0.4 mm/U
Overall lift 15 mm
Specimen stage 75 x 30 mm R/L with ceramic surface

Additional equipment (optional)

Objective Plan-Neofluar 20x/0.5; 40x/0.75
Digital camera (USB port) Lumenera 175c
Large bi-axial ceramic surface 102 x 105 mm R/L
Illumination gauge Blue
for measuring light intensity EUROIMMUN

EUROIMMUN, a specialist in the field of immunodiagnosics, now offers you the **EUROStar II fluorescence microscope**, developed in cooperation with a well-known manufacturer of optical equipment. This microscope is specifically tailored to the requirements of indirect immunofluorescence. Non-essential components have been removed and the labourious conventional mercury vapour lamp has been replaced by the stunningly simple EUROStar-Bluelight system.

Only a fraction of the light from a mercury lamp (HBO) can be used for fluorescence excitation – the majority of the energy is transformed into heat and dangerous ultraviolet radiation. With the EUROStar-Bluelight, engineers at EUROIMMUN AG have brought blue light-emitting diodes to fluorescence microscopy. Almost all the emitted light is suitable for the excitation of fluorescein.

Requiring only a tenth of the electrical power (low heat formation) the EUROStar-Bluelight provides the same amount of usable excitation light as a 50-watt HBO mercury vapour lamp.

The EUROStar-Bluelight does not emit any ultraviolet radiation and is explosion proof. We were therefore able to integrate this light source into the microscope housing, so there is no further need for the large

and often unreliable switching transformer. Also, there is no longer any need to keep records of lamp operation times – and these were in any case superfluous owing to the different operational life span of each individual HBO lamp.

With a frequently used fluorescence microscope, the mercury vapour lamp is often left switched on all day as, once switched off, it must be left to cool for 20 minutes before being switched on again, and several minutes are then needed after switching on until the full light output is achieved. Keeping a fluorescence microscope in constant operation in this way is costly in terms of energy use and seriously reduces the effective operating life of the lamp.

In contrast, the EUROStar-Bluelight provides instant full light output after being switched on, and can be switched off and on again in quick succession with no harmful effects. The LED has a life expectancy of at least 15,000 hours – in which time 150 mercury vapour lamps would have burned out, at great expense and requiring 150 beam adjustments to be made! When, at some point, the life of the LED ends, the operator is prompted in time by an acoustic alarm to arrange for an exchange of the light source.

EUROIMMUN has successfully overcome a serious weakness of the indirect immun-

ofluorescence: Depending on the design and age of the light source installed in the microscope which was used for evaluation, some antibodies had been evaluated as positive and, at other times, as negative or widely differing titers were determined for one and the same sample. The EUROStar II provides a defined light flux which is maintained at a constant level through the entire life span of the light source. The EUROStar II will help you to make an excellent impression during your next accreditation procedure.

Power for the EUROStar II microscope is supplied from a normal mains outlet (via the supplied adaptor, where applicable). Should a power loss occur, the microscope can even be powered for a short time from a medium-sized 12-volt battery!

Customers who require white transmitted light for conventional microscopy in addition to epifluorescence microscopy, can purchase the novel version **EUROStar II Plus**, which has an additional light source mounted on the condenser carrier to project a **cool, white LED light** directly under the object – the end of Köhler's illumination, one less hurdle in practical physics training to manage. Furthermore, the EUROStar II microscope has been designed in its standard to accept a digital camera, allowing immunofluorescence images to be captured and stored, or used for teaching and presentation purposes.